



# Stealth<sup>®</sup>

Atraumatic Occlusion System

## Advantages of the Stealth Clips

For occlusion where space is at a premium, the Stealth clip provides a 40% lower profile than our previous models.



### Reduced Spring Force

Three optimized spring forces to accommodate the procedure and the surgeon's preference.



### *latis*<sup>®</sup>/*latis* clips

Reorder No.	Description	Size
A1601	<i>latis/latis</i> clip	6mm (1/2 Force)
A1602	<i>latis/latis</i> clip	6mm (1/4 Force)
A1603	<i>latis/latis</i> clip	6mm (3/4 Force)

### Soft/Soft clips

Reorder No.	Description	Size
A1604	Soft/Soft clip	6mm (3/4 Force)





## Soft/Fibra® clips

Reorder No.	Description	Size
A1652	Soft/Fibra clip	6mm (1/2 Force)
A1654	Soft/Fibra clip	6mm (1/4 Force)
A1650	Soft/Fibra clip	6mm (3/4 Force)




## Stealth Surgical Clip Applicator

Reorder No.	Description	Handle Angle	Overall Length
A3200	Stealth Surgical Clip Applicator	13°	7.25"

## Advantages of the Stealth Occlusion System

Due to its superior traction, reduced occlusion force and minimized migration, the low profile Stealth system has become a leader in atraumatic occlusion.

## Comparison of the Stealth Clamp



## Ease of Installation

Secure and easy attachment



## Variety of Sizes

Stealth systems offer a wide variety of sizes and configurations of clamp inserts.

Size 1 - 30mm




Size 2 - 45mm




Size 3 - 60mm



Size 5 - 90mm



Sidebite - small  
50mm




Sidebite - medium  
75mm

*Actual sizes*

## Stealth Inserts

Stealth clamp inserts are available in *latis*, Fibra and Soft/Traction configurations to best fit your specific occlusion requirements.



### Soft Inserts

Our soft Stealth inserts are latex-free, as are all of our Stealth options.

### Soft/Soft Inserts

Reorder No.	Size
A0A11	30mm
A0A12	45mm
A0A13	60mm
A0A14	90mm



## ***latis* Inserts**

The unique *latis* mesh insert gently engages the adventitia reducing migration and compressive force.

Reorder No.	Size
-------------	------

### *Straight*

A0C01	30mm
A0C02	45mm
A0C03	60mm
A0C04	90mm

### *Sidebites*

A0C05	Small
A0C06	Medium



## **Soft/Traction Inserts**

Reorder No.	Size
-------------	------

### *Straight*

A0A01	30mm
A0A02	45mm
A0A03	60mm
A0A04	90mm



## **Fibra Inserts**

The Fibra inserts employ a woven layer of flexible, soft, finger-like fibers that mesh with the adventitia to provide pressure-independent traction.

## **Fibra/Fibra Inserts**

Reorder No.	Size
-------------	------

### *Straight*

A0G11	30mm
A0G12	45mm
A0G13	60mm
A0G14	90mm

### *Sidebites*

A0G05	Small
A0G06	Medium



## Soft/Fibra Inserts

Reorder No.	Size
-------------	------

*Straight*

A0G01	30mm
-------	------

A0G02	45mm
-------	------

A0G03	60mm
-------	------

A0G04	90mm
-------	------

## Stealth Clamps

Atraumatic Stealth clamps offer an unmatched level of versatility.

### Consistent, Gentle Occlusion


Parallel jaw alignment provides consistent, gentle occlusion

### Low Profile

Optimal occlusion when exposure is limited

### Ultimate Insert Security

Advanced attachment provides ultimate security of insert on clamp





### Ensured Durability

Stainless Steel Finish

### Enhanced Stability

Exceptionally crafted fine-tuned box lock improves stability

## Preferred Configurations



### Bahnsen

Stealth system's best-selling model



### Wylie

Deep reach with unique bend



### Sidebite

The market's only padded sidebite

## Stealth Clamp Configurations



Reorder no.	Size	Handle angle	Overall length	Jaw angle	Jaw length	Applications
A3201	1	15°	6.5"	Straight	30mm	Femoral and Iliac
A3202	1	25°	6"	60°	30mm	Femoral and Iliac

## Stealth Clamp Configurations



Reorder no.	Size	Handle angle	Overall length	Jaw angle	Jaw length	Applications
A3203	2	15°	7"	Straight	45mm	Femoral and Iliac
A3204	2	15°	7"	34°	45mm	Femoral and Iliac
A3205	2	Straight	6.5"	65°	45mm	Femoral and Iliac
A3206	2	15°	6.5"	65°	45mm	Femoral and Iliac
A3207	3	10°	11.5"	Straight	60mm	Iliac and Abdominal Aortic
A3208	3	15°	10"	60°	60mm	Iliac and Abdominal Aortic
A3209	3	35°	10"	60°	60mm	Iliac and Abdominal Aortic
A3210	3	35°	10"	Bahnson	60mm	Iliac and Abdominal Aortic
A3211	3	60°	10.5"	Wylie	60mm	Iliac, Abdominal Aortic, Femoral and Renal
A3212	5	15°	12"	Straight	90mm	Aortic Cross Clamp
A3213	5	35°	11"	60°	90mm	Aortic Cross Clamp
A3214	5	35°	10.5"	Bahnson	90mm	Aortic Cross Clamp
A3215	Small	25°	7"	Sidebite	50mm	Aortic and Vena Cava
A3216	Medium	45°	8"	Sidebite	75mm	Aortic and Vena Cava
A3219	3	35°	8.5"	Bahnson	60mm	Aortic Cross Clamp
A3220	5	35°	9"	Bahnson	90mm	Aortic Cross Clamp
A3221	2	8.5°	7.5"	Debakey	45mm	Femoral and Iliac

Visit [www.appliedmedical.com/stealth](http://www.appliedmedical.com/stealth) for more information

Devices listed may not be approved in all markets.  
Please contact your Field Implementation team member for more information on availability.

© 2019 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union.  
213951-EN-USA-B

